

FIXATION SYSTEMS



MB1-17-001-45

MF-height adjustable headholder, including 3 gel pads

- ✓ Universally applicable with 3 gel pads for safe positioning of the head
- ✓ Particularly suitable for newborns, children and adolescents, individually adjustable from a head diameter of 10cm
- ✓ Can be used as head support or together with skull clamp
- ✓ MB1-17-001-46: Gel pad for MB1-17-001-45, please order 3 pcs in case of replacement
- ✓ To be used with MF-Swivel adapter Tripod MB1-17-001-30

MB1-17-001-45

MF-Höhenverstellbare Kopfstütze, inklusive 3 Gelpads

- ✓ Universell einsetzbar mit 3 Gelpads zur sicheren Lagerung des Kopfes
- ✓ Besonders geeignet für Neugeborene, Kinder und Jugendliche, individuell einstellbar ab einem Kopfdurchmesser von 10cm
- ✓ Als Kopflagerung oder zusammen mit Schädelklemme einsetzbar
- ✓ MB1-17-001-46: Gelpad für MB1-17-001-45, bitte 3 Stück bestellen als Ersatz
- ✓ Zu verwenden mit MF-Drehadapter MB1-17-001-30



Mehr zu unserem MF-Kopfhalterungssystem finden Sie hier.
Scannen oder klicken



For more information about our MF-head-rest Cranial fixation system
scan or click